



KCNIP3 Antibody, HRP conjugated

Product Code	CSB-PA012045LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2W7
Immunogen	Recombinant Human Calsenilin protein (1-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calsenilin (A-type potassium channel modulatory protein 3) (DRE-antagonist modulator) (DREAM) (Kv channel-interacting protein 3) (KCHIP3), KCNIP3, CSEN DREAM KCHIP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	KCNIP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.