



KCNIP3 Antibody, Biotin conjugated

Product Code	CSB-PA012045LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2W7
Immunogen	Recombinant Human Calsenilin protein (1-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calsenilin (A-type potassium channel modulatory protein 3) (DRE-antagonist modulator) (DREAM) (Kv channel-interacting protein 3) (KCHIP3), KCNIP3, CSEN DREAM KCHIP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	KCNIP3