



# KCNK1 Antibody, FITC conjugated

Product Code	CSB-PA012062LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00180
Immunogen	Recombinant Human Potassium channel subfamily K member 1 protein (268-336AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel subfamily K member 1 (Inward rectifying potassium channel protein TWIK-1) (Potassium channel K2P1) (Potassium channel KCNO1), KCNK1, HOHO1 KCNO1 TWIK1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNK1