



KIF20A Antibody, HRP conjugated

Product Code	CSB-PA012321LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95235
Immunogen	Recombinant Human Kinesin-like protein KIF20A protein (161-294AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF20A (GG10_2) (Mitotic kinesin-like protein 2) (MKlp2) (Rab6-interacting kinesin-like protein) (Rabkinesin-6), KIF20A, MKLP2, RAB6KIFL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIF20A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.