



# KLF7 Antibody, HRP conjugated

Product Code	CSB-PA012397LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75840
Immunogen	Recombinant Human Krueppel-like factor 7 protein (114-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Krueppel-like factor 7 (Ubiquitous krueppel-like factor), KLF7, UKLF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KLF7