



KLK8 Antibody, Biotin conjugated

Product Code	CSB-PA012459LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60259
Immunogen	Recombinant Human Kallikrein-8 protein (87-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-8 (hK8) (EC 3.4.21.118) (Neuropsin) (NP) (Ovasin) (Serine protease 19) (Serine protease TADG-14) (Tumor-associated differentially expressed gene 14 protein), KLK8, NRPN PRSS19 TADG14
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KLK8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.