



MEOX2 Antibody, FITC conjugated

Product Code	CSB-PA013696LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50222
Immunogen	Recombinant Human Homeobox protein MOX-2 protein (94-193AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein MOX-2 (Growth arrest-specific homeobox) (Mesenchyme homeobox 2), MEOX2, GAX MOX2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MEOX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.