



MNAT1 Antibody, HRP conjugated

Product Code	CSB-PA014686LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51948
Immunogen	Recombinant Human CDK-activating kinase assembly factor MAT1 protein (1-309AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CDK-activating kinase assembly factor MAT1 (CDK7/cyclin-H assembly factor) (Cyclin-G1-interacting protein) (Menage a trois) (RING finger protein 66) (RING finger protein MAT1) (p35) (p36), MNAT1, CAP35 MAT1 RNF66
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MNAT1