



NCR1 Antibody, Biotin conjugated

Product Code	CSB-PA015549LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76036
Immunogen	Recombinant Human Natural cytotoxicity triggering receptor 1 protein (121-230AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural cytotoxicity triggering receptor 1 (Lymphocyte antigen 94 homolog) (NK cell-activating receptor) (Natural killer cell p46-related protein) (NK-p46) (NKp46) (hNKp46) (CD antigen CD335), NCR1, LY94
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NCR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.