



NQO2 Antibody, FITC conjugated

Product Code	CSB-PA016040LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16083
Immunogen	Recombinant Human Ribosyldihyronicotinamide dehydrogenase [quinone] protein (1-231AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NQO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.