



NQO2 Antibody, Biotin conjugated

Product Code	CSB-PA016040LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16083
Immunogen	Recombinant Human Ribosyldihyronicotinamide dehydrogenase [quinone] protein (1-231AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosyldihyronicotinamide dehydrogenase [quinone] (EC 1.10.5.1) (NRH dehydrogenase [quinone] 2) (NRH:quinone oxidoreductase 2) (Quinone reductase 2) (QR2), NQO2, NMOR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NQO2