




RPS27A Antibody

Product Code	CSB-PA01639A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62979
Immunogen	Recombinant Human Ubiquitin-40S ribosomal protein S27a protein (16-56AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-40S ribosomal protein S27a (Ubiquitin carboxyl extension protein 80) [Cleaved into: Ubiquitin; 40S ribosomal protein S27a (Small ribosomal subunit protein eS31)], RPS27A, UBA80 UBCEP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RPS27A

Image



Immunofluorescent analysis of MCF-7 cells using CSB-PA01639A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.