



RPS27A Antibody, Biotin conjugated

Product Code	CSB-PA01639D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62979
Immunogen	Recombinant Human Ubiquitin-40S ribosomal protein S27a protein (16-56AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-40S ribosomal protein S27a (Ubiquitin carboxyl extension protein 80) [Cleaved into: Ubiquitin; 40S ribosomal protein S27a (Small ribosomal subunit protein eS31)], RPS27A, UBA80 UBCEP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RPS27A