




PPIC Antibody

Product Code	CSB-PA018473LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45877
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase C protein (26-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase C (PPIase C) (EC 5.2.1.8) (Cyclophilin C) (Rotamase C), PPIC, CYPC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PPIC

Image



Immunofluorescent analysis of HeLa cells using CSB-PA018473LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)