



PPIC Antibody, Biotin conjugated

Product Code	CSB-PA018473LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45877
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase C protein (26-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase C (PPIase C) (EC 5.2.1.8) (Cyclophilin C) (Rotamase C), PPIC, CYPC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PPIC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.