



PRDX6 Antibody, Biotin conjugated

Product Code	CSB-PA018659LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30041
Immunogen	Recombinant Human Peroxiredoxin-6 protein (87-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxiredoxin-6 (EC 1.11.1.15) (1-Cys peroxiredoxin) (1-Cys PRX) (24 kDa protein) (Acidic calcium-independent phospholipase A2) (aiPLA2) (EC 3.1.1.4) (Antioxidant protein 2) (Liver 2D page spot 40) (Non-selenium glutathione peroxidase) (NSGPx) (Red blood cells page spot 12), PRDX6, AOP2 KIAA0106
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRDX6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.