



RORA Antibody, HRP conjugated

Product Code	CSB-PA020069LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35398
Immunogen	Recombinant Human Nuclear receptor ROR-alpha protein (110-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor ROR-alpha (Nuclear receptor RZR-alpha) (Nuclear receptor subfamily 1 group F member 1) (RAR-related orphan receptor A) (Retinoid-related orphan receptor-alpha), RORA, NR1F1 RZRA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RORA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.