



RORA Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA020069LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P35398 |
| Immunogen | Recombinant Human Nuclear receptor ROR-alpha protein (110-214AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nuclear receptor ROR-alpha (Nuclear receptor RZR-alpha) (Nuclear receptor subfamily 1 group F member 1) (RAR-related orphan receptor A) (Retinoid-related orphan receptor-alpha), RORA, NR1F1 RZRA |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | RORA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |