



SERINC1 Antibody, FITC conjugated

Product Code	CSB-PA021046LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRX5
Immunogen	Recombinant Human Serine incorporator 1 protein (331-387AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine incorporator 1 (Tumor differentially expressed protein 1-like) (Tumor differentially expressed protein 2), SERINC1, KIAA1253 TDE1L TDE2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SERINC1