



SERINC1 Antibody, Biotin conjugated

Product Code	CSB-PA021046LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRX5
Immunogen	Recombinant Human Serine incorporator 1 protein (331-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine incorporator 1 (Tumor differentially expressed protein 1-like) (Tumor differentially expressed protein 2), SERINC1, KIAA1253 TDE1L TDE2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SERINC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.