




SPINT2 Antibody

Product Code	CSB-PA022585LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43291
Immunogen	Recombinant Human Kunitz-type protease inhibitor 2 protein (71-194AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kunitz-type protease inhibitor 2 (Hepatocyte growth factor activator inhibitor type 2) (HAI-2) (Placental bikunin), SPINT2, HAI2 KOP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SPINT2

Image



Western Blot

Positive WB detected in: Mouse lung tissue


All lanes: SPINT2 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 29, 22 kDa

Observed band size: 29 kDa



Immunofluorescent analysis of A549 cells using CSB-PA022585LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.