



SPINT2 Antibody, Biotin conjugated

Product Code	CSB-PA022585LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43291
Immunogen	Recombinant Human Kunitz-type protease inhibitor 2 protein (71-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kunitz-type protease inhibitor 2 (Hepatocyte growth factor activator inhibitor type 2) (HAI-2) (Placental bikunin), SPINT2, HAI2 KOP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SPINT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.