



TAT Antibody, Biotin conjugated

Product Code	CSB-PA023175LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17735
Immunogen	Recombinant Human Tyrosine aminotransferase protein (25-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine aminotransferase (TAT) (EC 2.6.1.5) (L-tyrosine:2-oxoglutarate aminotransferase), TAT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TAT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.