



TBL1X Antibody, Biotin conjugated

Product Code	CSB-PA023234LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60907
Immunogen	Recombinant Human F-box-like/WD repeat-containing protein TBL1X protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box-like/WD repeat-containing protein TBL1X (SMAP55) (Transducin beta-like protein 1X) (Transducin-beta-like protein 1, X-linked), TBL1X, TBL1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TBL1X
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.