



TOM1 Antibody, FITC conjugated

Product Code	CSB-PA024043LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60784
Immunogen	Recombinant Human Target of Myb protein 1 protein (306-487AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Target of Myb protein 1, TOM1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TOM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.