




UBE2L6 Antibody

Product Code	CSB-PA025466LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14933
Immunogen	Recombinant Human Ubiquitin/ISG15-conjugating enzyme E2 L6 protein (13-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin/ISG15-conjugating enzyme E2 L6 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme L6) (Retinoic acid-induced gene B protein) (RIG-B) (Ubch8) (Ubiquitin carrier protein L6) (Ubiquitin-protein ligase L6), UBE2L6, UBCH8
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2L6

Image



Immunofluorescent analysis of HeLa cells using CSB-PA025466LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.