



USO1 Antibody, Biotin conjugated

Product Code	CSB-PA025694LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60763
Immunogen	Recombinant Human General vesicular transport factor p115 protein (591-733AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	General vesicular transport factor p115 (Protein USO1 homolog) (Transcytosis-associated protein) (TAP) (Vesicle-docking protein), USO1, VDP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	USO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.