




TNFRSF11B Antibody

Product Code	CSB-PA08229A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00300
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 11B protein (22-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 11B (Osteoclastogenesis inhibitory factor) (Osteoprotegerin), TNFRSF11B, OCIF OPG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TNFRSF11B

Image



Immunohistochemistry of paraffin-embedded human prostate tissue using CSB-PA08229A0Rb at dilution of 1:100