



TNFRSF11B Antibody, Biotin conjugated

Product Code	CSB-PA08229D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00300
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 11B protein (22-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 11B (Osteoclastogenesis inhibitory factor) (Osteoprotegerin), TNFRSF11B, OCIF OPG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	TNFRSF11B