



CLDN3 Antibody, Biotin conjugated

Product Code	CSB-PA13119D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15551
Immunogen	Recombinant Human Claudin-3 protein (181-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Claudin-3 (Clostridium perfringens enterotoxin receptor 2) (CPE-R 2) (CPE-receptor 2) (Rat ventral prostate.1 protein homolog) (hRVP1), CLDN3, C7orf1 CPETR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CLDN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.