





CRAT Antibody

Product Code	CSB-PA14169A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43155
Immunogen	Recombinant Human Carnitine O-acetyltransferase protein (181-346AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carnitine O-acetyltransferase (Carnitine acetylase) (EC 2.3.1.7) (Carnitine acetyltransferase) (CAT) (CrAT), CRAT, CAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CRAT

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA14169A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA14169A0Rb at dilution of 1:100