



MSH2 Antibody, FITC conjugated

Product Code	CSB-PA16889C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43246
Immunogen	Recombinant Human DNA mismatch repair protein Msh2 protein (2-254AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA mismatch repair protein Msh2 (hMSH2) (MutS protein homolog 2), MSH2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MSH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.