




omcB Antibody

Product Code	CSB-PA315528HA01DSB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CC04
Immunogen	Recombinant Chlamydia trachomatis Large cysteine-rich periplasmic protein OmcB protein (15-243AA)
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Large cysteine-rich periplasmic protein OmcB (Large-CRP) (60 kDa cysteine-rich OMP) (60 kDa outer membrane protein) (Cysteine-rich outer membrane protein) (CRP), omcB, omp2 omp2B
Immunogen Species	Chlamydia trachomatis
Research Area	Others
Target Names	omcB

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: omcB antibody at 3.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.