



yebF Antibody, Biotin conjugated

Product Code	CSB-PA329240HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33219
Immunogen	Recombinant Escherichia coli Protein YebF protein (22-118AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein YebF, yebF
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	yebF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.