



yebF Antibody, Biotin conjugated

Product Code	CSB-PA329240HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33219
Immunogen	Recombinant Escherichia coli Protein YebF protein (22-118AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein YebF, yebF
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	yebF