




ftsZ Antibody

Product Code	CSB-PA359270HA01EGX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A9A7
Immunogen	Recombinant Escherichia coli O6:H1 Cell division protein FtsZ protein (1-383AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O6:H1
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division protein FtsZ, ftsZ
Immunogen Species	Escherichia coli O6:H1
Research Area	Others
Target Names	ftsZ

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ftsZ antibody at 2.7μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 58 kDa

Observed band size: 58 kDa