



ftsZ Antibody, Biotin conjugated

Product Code	CSB-PA359270HD01EGX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A9A7
Immunogen	Recombinant Escherichia coli O6:H1 Cell division protein FtsZ protein (1-383AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O6:H1
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division protein FtsZ, ftsZ
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Research Area	Others
Target Names	ftsZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.