



fbpC Antibody, FITC conjugated

Product Code	CSB-PA363614LC01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WQN8
Immunogen	Recombinant Mycobacterium tuberculosis Diacylglycerol acyltransferase/mycolyltransferase Ag85C protein (46-340AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diacylglycerol acyltransferase/mycolyltransferase Ag85C (DGAT) (EC 2.3.1.122) (EC 2.3.1.20) (Acyl-CoA:diacylglycerol acyltransferase) (Antigen 85 complex C) (85C) (Ag85C) (Fibronectin-binding protein C) (Fbps C), fbpC, mpt45
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Research Area	Others
Target Names	fbpC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.