



MIA Antibody, HRP conjugated

Product Code	CSB-PA613694LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16674
Immunogen	Recombinant Human Melanoma-derived growth regulatory protein (87-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-derived growth regulatory protein (Melanoma inhibitory activity protein), MIA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.