



MIA Antibody, HRP conjugated

Product Code	CSB-PA613694LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16674
Immunogen	Recombinant Human Melanoma-derived growth regulatory protein (87-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Melanoma-derived growth regulatory protein (Melanoma inhibitory activity protein), MIA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MIA