



DPYSL2 Antibody, Biotin conjugated

Product Code	CSB-PA614975YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16555
Immunogen	Recombinant Human Dihydropyrimidinase-related protein 2 protein (1-572AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydropyrimidinase-related protein 2 (DRP-2) (Collapsin response mediator protein 2) (CRMP-2) (N2A3) (Unc-33-like phosphoprotein 2) (ULIP-2), DPYSL2, CRMP2 ULIP2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DPYSL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.