



NDUFA5 Antibody, HRP conjugated

Product Code	CSB-PA614988LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16718
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 5 protein (38-114AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 5 (Complex I subunit B13) (Complex I-13kD-B) (CI-13kD-B) (NADH-ubiquinone oxidoreductase 13 kDa-B subunit), NDUFA5
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NDUFA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.