




CDH18 Antibody

Product Code	CSB-PA615709LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13634
Immunogen	Recombinant Human Cadherin-18 protein (66-251AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin-18 (Cadherin-14), CDH18, CDH14
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CDH18

Image



Western Blot

Positive WB detected in: Mouse kidney tissue


All lanes: CDH18 antibody at 2.4μg/ml

Secondary


Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 89, 64 kDa


Observed band size: 89 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA615709LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA615709LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA615709LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)