



PKN2 Antibody, HRP conjugated

Product Code	CSB-PA620966LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16513
Immunogen	Recombinant Human Serine/threonine-protein kinase N2 protein (925-969AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase N2 (EC 2.7.11.13) (PKN gamma) (Protein kinase C-like 2) (Protein-kinase C-related kinase 2), PKN2, PRK2 PRKCL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PKN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.