



TRIP10 Antibody, Biotin conjugated

Product Code	CSB-PA624029LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15642
Immunogen	Recombinant Human Cdc42-interacting protein 4 protein (276-486AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cdc42-interacting protein 4 (Protein Felic) (Salt tolerant protein) (hSTP) (Thyroid receptor-interacting protein 10) (TR-interacting protein 10) (TRIP-10), TRIP10, CIP4 STOT STP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIP10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.