



TFDP3 Antibody, Biotin conjugated

Product Code	CSB-PA702134LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5H9I0
Immunogen	Recombinant Human Transcription factor Dp family member 3 protein (61-196AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription factor Dp family member 3 (Cancer/testis antigen 30) (CT30) (Hepatocellular carcinoma-associated antigen 661), TFDP3, DP4 HCA661
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TFDP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.