



MICU2 Antibody, Biotin conjugated

Product Code	CSB-PA818265LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYU8
Immunogen	Recombinant Human Calcium uptake protein 2, mitochondrial protein (26-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium uptake protein 2, mitochondrial (EF-hand domain-containing family member A1), MICU2, EFHA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MICU2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.