



MICU2 Antibody, Biotin conjugated

Product Code	CSB-PA818265LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IYU8
Immunogen	Recombinant Human Calcium uptake protein 2, mitochondrial protein (26-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium uptake protein 2, mitochondrial (EF-hand domain-containing family member A1), MICU2, EFHA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MICU2