



SLC20A1 Antibody, HRP conjugated

Product Code	CSB-PA823882LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUM9
Immunogen	Recombinant Human Sodium-dependent phosphate transporter 1 protein (263-458AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent phosphate transporter 1 (Gibbon ape leukemia virus receptor 1) (GLVR-1) (Leukemia virus receptor 1 homolog) (Phosphate transporter 1) (PiT-1) (Solute carrier family 20 member 1), SLC20A1, GLVR1 PIT1
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	SLC20A1