



SLC20A1 Antibody, Biotin conjugated

Product Code	CSB-PA823882LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUM9
Immunogen	Recombinant Human Sodium-dependent phosphate transporter 1 protein (263-458AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent phosphate transporter 1 (Gibbon ape leukemia virus receptor 1) (GLVR-1) (Leukemia virus receptor 1 homolog) (Phosphate transporter 1) (PiT-1) (Solute carrier family 20 member 1), SLC20A1, GLVR1 PIT1
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	SLC20A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.